

Cap-ONE CO₂ Sensor TG-920P



A Whole
New Way of
Looking at
Capnography
Monitoring

- Mainstream technology CO₂ monitoring for both intubated and non-intubated patients
- Smart modular cable technology is compatible with Nihon Kohden bedside monitors and stand-alone Cap-STAT monitors
- Simple operation with no warm-up time or calibration needed

Cap-ONE CO₂ Sensor TG-920P

If you're prescribing CO₂ monitoring, shouldn't it work under all conditions?

Whether your patients are oral or nasal breathers, in need of short or long-term monitoring or under high-humidity conditions, they require and deserve reliable capnography monitoring.

Now you can achieve the results you need with Cap-ONE, The world's 1st mainstream CO₂ sensor specifically designed for intubated or non-intubated patients.

Using advanced miniaturization and sensor technology, we have substantially reduced the size of our traditional mainstream sensor.


This new sensor is attached to a disposable oral and nasal adaptor and placed directly at the point of expiration.

Therefore, you can achieve the same level of quality and reliability found in traditional mainstream CO₂ technology monitoring and apply these benefits to intubated or non-intubated patients without any of the hassle and cross-contamination concerns found in traditional sidestream technology.

Clinical research has demonstrated that Cap-ONE proves to be more effective in capnography monitoring for non-intubated patients than traditional sidestream technology, as shown below.



Cap-ONE delivers safe and accurate capnography monitoring with virtually no delay in data sampling and minimal discomfort to your patients.



Patient starts breathing partially with mouth.