



NK-HiQ™ CNS Remote Client vs Life Scope® Central Station

MONITORING


FEATURES	Central Station	CNS Remote Client *
4, 6, 8, 10, 12 and 16 Patient Views	Yes	Yes
Supports Dual Screen Configuration	Yes	No
Primary Bed Monitor (local file storage)	Yes	No
Touch Screen Operation	Yes	No
Audible and Visual Alarms	Yes	Yes
Use Alarms for Primary Alarm Notification	Yes	No
Silence Alarms	Yes	Yes
All Beds Customization “Drag & Drop”	Yes	Yes
Change Alarm Limits	Yes	Yes

REVIEW

FEATURES	Central Station	CNS Remote Client *
View Waveforms	Yes	Yes
View Vital Signs	Yes	Yes
View Trend Data	Yes	Yes
View Full Disclosure	Yes	Yes
View Arrhythmia Recall	Yes	Yes
View Event List	Yes	Yes
View 12-Lead ECGs	Yes	Yes
View All Leads of ECG	Yes	Yes
View ST Recall	Yes	Yes
View Hemodynamics	Yes	Yes

FEATURES

FEATURES	Central Station	CNS Remote Client *
Change Waveform Size	Yes	Yes
Send ECG Strips to EMR Support	Yes	Yes
Print All	Yes	Yes
Strip Recording	Yes	Yes
Admit/Discharge/Transfer	Yes	Yes
Data Archival/Readmit	Yes	No
Change Device (when bed is local stored)	Yes	N/A
Change Remote Filing Bed (when not local stored)	No	Yes
Change Patient Demographic Information	Yes	Yes
Change Transmitter Channel	Yes	Yes
Pause Monitoring	Yes	Yes
Adjust Alarm Volume	Yes	Yes
Set Bed Color	Yes	Yes



* (Used for secondary monitoring)

NK-HiQ is a trademark of Nihon Kohden Corporation.
Life Scope and Different Thinking for Better Healthcare are registered trademarks of Nihon Kohden.

Different Thinking for Better Healthcare.®

