



# NK-HiQ™ CNS Remote Client vs Remote Station

## MONITORING


FEATURES	CNS Remote Client *	Remote Station *
<b>Audible and Visual Alarms</b>	Yes	Yes
<b>Use Alarms for Secondary Alarm Notification</b>	Yes	Yes
<b>Silence Alarms</b>	Yes	Yes
<b>Change Alarm Limits</b>	Yes	Yes

## REVIEW

FEATURES	CNS Remote Client *	Remote Station *
<b>4, 6, 8, 10, 12 and 16 Patient Views</b>	Yes	Yes
<b>Maximum number of beds per display</b>	16	16
<b>Primary Bed Monitor (local file storage)</b>	No	No
<b>Touch Screen Operation</b>	No	No
<b>All Beds Customization “Drag &amp; Drop”</b>	Yes	No
<b>Review Period</b>	120 hours	72 hours
<b>View Vital Signs</b>	Yes	Yes
<b>View Trend Data</b>	Yes	Yes
<b>View Full Disclosure</b>	Yes Up to 16 Waves	Yes Up to 6 Waves
<b>View Arrhythmia Recall</b>	Yes	Yes
<b>View Event List (Alarm History)</b>	Yes	No
<b>View 12 Lead ECGs</b>	Yes	Yes
<b>View Alarm Events Window</b>	Yes	N/A
<b>View ST Recall</b>	Yes	Yes
<b>View Hemodynamics</b>	Yes	Yes
<b>View SpO<sub>2</sub> Trend</b>	Yes	N/A

## FEATURES

FEATURES	CNS Remote Client *	Remote Station *
<b>Send ECG Strips to EMR Support</b>	Yes	No
<b>Print All</b>	Yes	No
<b>Strip Recording</b>	Yes	Yes
<b>Admit/Discharge/Transfer</b>	Yes	Yes
<b>Change Patient Demographic Information</b>	Yes	Yes
<b>Pause Monitoring</b>	Yes	Yes
<b>Data Archival/Readmit</b>	No	No
<b>End User Bed Assignment</b>	Yes	Yes
<b>Change Transmitter Channel</b>	Yes	Yes
<b>Change Device</b>	No	No
<b>Adjust Alarm Volume</b>	Yes	Yes
<b>Set Bed Color</b>	Yes	Yes



\* (Used for secondary monitoring)

NK-HiQ is a trademark of Nihon Kohden Corporation.  
Different Thinking for Better Healthcare is a registered trademark of Nihon Kohden.

Different Thinking for  
Better Healthcare.®

