

VITRAC[®]

REMOTE VIEWING APP



Connect Anywhere, Anytime



Data at Your Fingertips

- Provides a secure method for viewing a range of Nihon Kohden-generated patient data that empowers clinicians to make timely decisions and supports quality of care
- Offers a continuous and simultaneous display of data for up to 12 patients on iPad[®] and Microsoft Windows[®] devices and up to eight patients on iPhone[®], Google Pixel[™] 1, Samsung Galaxy[®] Note8, and Zebra[®] TC51 Android[™] devices
- Offers a robust, easy-to-use interface, which allows users to see current waveforms, vital signs, stored data, and much more
- Compatible with Apple[®] mobile iOS devices, Android devices, and Microsoft Windows[®] devices

Interactive Monitoring Data

- Graphical and tabular trends
- Hemodynamic calculations
- Arrhythmia events
- Minute-to-minute ST templates
- Multi-parameter full disclosure waveforms
- Interpretative 12-lead ECGs

Secure Viewing

- Secured access for local or remote users
- Near real-time waveforms and numerics with user-selectable settings for customized views
- Allows data to be viewed in near real-time within the hospital network or remotely via a VPN connection
- Leverages NK-HiQ[™] Enterprise Gateway multi-level security controls, which restrict access to patient data based on user, user groups, and/or mobile client

VITRAC® REMOTE VIEWING APP


Compatible with Multiple Devices

- Apple iPad Pro 12.9", running iOS 11 or higher
- Apple iPad Mini 4, running iOS 11 or higher
- Apple iPhone 7, 8, X, XS, running iOS 11 or higher
- Google Pixel 1, Samsung Galaxy Note 8, and Zebra TC51 running Android 4.1 (Jelly Bean) and higher, up-to Android 8.1 (Oreo)
- Microsoft OS devices running Windows 10

System Overview

The ViTrac remote viewing application is part of the NK-HiQ Enterprise Gateway platform. Designed for secondary viewing, the application includes all required software to provide clinicians with secure remote access to patient data.

As the backbone of the Nihon Kohden Patient Monitoring System, the NK-HiQ Enterprise Gateway platform consolidates multiple applications into a single powerful server — saving time, space, and money.



For more information, please contact us at 1-800-325-0283 or visit us.nihonkohden.com



Exceptional Solutions,
Simply Delivered.™

ViTrac is a registered trademark of NKUS Lab, a subsidiary of Nihon Kohden Corporation. NK-HiQ is a trademark of Nihon Kohden Corporation. Apple, iPad, iPhone and iPod are registered trademarks of Apple Inc. Samsung Galaxy is a trademark of Samsung. Zebra is a trademark of ZIH Corp., registered in many jurisdictions worldwide. Google Pixel and Android are trademarks of Google LLC. Windows is a trademark of Microsoft Corporation. Exceptional Solutions, Simply Delivered is a trademark of Nihon Kohden.