

# CENTRAL STATION vs NGNK Remote Monitoring Station

MONITORING	CNS-6801	*NGNK-6803/04
<b>Models</b>	Central Station	Remote Monitoring Station
<b>4,6,8,10,12 and 16 Patient Views</b>	Yes	Yes
<b>Supports Dual Screen Configuration</b>	Yes	No
<b>Can be used as primary file storage (in addition to bedside)</b>	Yes	No
<b>Touch Screen Operation</b>	Yes	Optional
<b>Audible and Visual Alarms</b>	Yes	Yes
<b>Can be used as Primary Alarm Notification</b>	Yes	No
<b>Silence Alarms</b>	Yes	Yes
<b>All Beds Customization "Drag &amp; Drop"</b>	Yes	Yes
<b>Change Alarm Limits</b>	Yes	Yes

## REVIEW

Models	Central Station	* Remote Monitoring Station
<b>View Waveforms</b>	Yes	Yes
<b>View Vital Signs</b>	Yes	Yes
<b>View Trend Data</b>	Yes	Yes
<b>View Full Disclosure</b>	Yes	Yes
<b>View Arrhythmia Recall</b>	Yes	Yes
<b>View Event List</b>	Yes	Yes
<b>View 12-Lead ECGs</b>	Yes	Yes
<b>View All Leads of ECG</b>	Yes	Yes
<b>View ST Recall</b>	Yes	Yes
<b>View Hemodynamics</b>	Yes	Yes

## FEATURES

Models	Central Station	* Remote Monitoring Station
<b>Change Waveform Size</b>	Yes	Yes
<b>Send ECG Strips to EMR Support</b>	Yes	Yes
<b>Print All</b>	Yes	Yes
<b>Strip Recording</b>	Yes	Yes
<b>Admit/Discharge/Transfer</b>	Yes	Yes
<b>Data Archival/Readmit</b>	Yes	No
<b>Change Device (when bed is local stored)</b>	Yes	N/A
<b>Change Remote Filing Bed (when not local stored)</b>	No	Yes
<b>Change Patient Demographic Information</b>	Yes	Yes
<b>Change Transmitter Channel</b>	Yes	Yes
<b>Pause Monitoring</b>	Yes	Yes
<b>Adjust Alarm Volume</b>	Yes	Yes
<b>Set Bed Color</b>	Yes	Yes
<b>Secure Communication – Encrypted Data Channels</b>	No	Yes

\* Used for secondary monitoring only.