

CENTRAL STATION VS NGNK REMOTE MONITORING STATION

| | CNS-2101 | *NGNK-6803/04 |
|---|------------------|----------------------------|
| MONITORING | | |
| Models | Central Station | Remote Monitoring Station |
| 2, 4, 6, 8, 10, 12 and 16 Patient Views | Yes | Yes |
| Supports Dual Screen Configuration | Yes ¹ | No |
| Can be used as primary file storage (in addition to bedside) | Yes | No |
| Touch Screen Operation | Yes | Yes |
| Audible and Visual Alarms | Yes | 6803 – Yes / 6804 – No |
| Can be used as Primary Alarm Notification | Yes | No |
| Silence Alarms | Yes | Yes |
| All Beds Customization “Drag & Drop” | Yes | Yes |
| Change Alarm Limits | Yes | Yes |
| REVIEW | | |
| Models | Central Station | *Remote Monitoring Station |
| View Waveforms | Yes | Yes |
| View Vital Signs | Yes | Yes |
| View Trend Data | Yes | Yes |
| View Full Disclosure | Yes | Yes |
| View Arrhythmia Recall | Yes | Yes |
| View Event List | Yes | Yes |
| View 12-Lead ECGs | Yes | Yes |
| View All Leads of ECG | Yes | Yes |
| View ST Recall | Yes | Yes |
| View Hemodynamics | Yes | Yes |
| FEATURES | | |
| Models | Central Station | *Remote Monitoring Station |
| Change Waveform Size | Yes | Yes |
| Send ECG Strips to EMR Support | Yes | Yes |
| Print All | Yes | Yes |
| Strip Recording | Optional | No |
| Admit/Discharge/Transfer | Yes | Yes |
| Data Archival/Readmit | Yes | No |
| Change Device (when bed is local stored) | Yes | N/A |
| Change Remote Filing Bed (when not local stored) | No | Yes |
| Change Patient Demographic Information | Yes | Yes |
| Change Transmitter Channel | Yes | Yes |
| Pause Monitoring | Yes | Yes |
| Adjust Alarm Volume | Yes | Yes |
| Set Bed Color | Yes | Yes |
| Secure Communication – Encrypted Data Channels | Yes | Yes |

* Used for secondary monitoring only.

¹ The use of the expansion capability removes the ability to drive external displays such as Non-interactives or KVM solutions.

CAUTION: Federal (United States) Law restricts this device to sale by or on the order of a physician. See Operators Manual for full prescribing information, including indications for use, contraindications, warnings, precautions and adverse events.