

Central Nurse Station

- Continuous monitoring of system alerts to identify potential issues and threats
- Automatic detection, notifications and alerts that reduce and may help prevent system failures
- Advanced protection against cybersecurity threats
- Proactive notifications that help optimize performance and system uptime

Enterprise Gateway & Host 1000 Server

- Real-time system monitoring to maintain patient data flow
- Provide advanced security, analysis, and notifications of potential threats
- Multi-layered protection against cybersecurity threats like malware and ransomware

NGNK Remote Client

- Continuous system health monitoring to prevent disruptions
- Performance optimization for reliable remote access

Network Switches

- Automatic detection of hardware issues on switches provided by Nihon Kohden
- Network performance optimization for reliable data transmission

NK NetworkCare™ Server

- Continuous monitoring to ensure system stability and reliability
- Performance optimization for seamless operation
- Malware and ransomware protection

Bedside Monitors


- Optional protection against cybersecurity threats



CNS-2101 Central Nurse Station



Enterprise Gateway Server



Life Scope G7 Bedside Monitor

NIHON KOHDEN AMERICA
15353 Barranca Pkwy, Irvine, CA 92618
800.325.0283 us.nihonkohden.com

Important Safety Information:

Caution: Federal (United States) law restricts this device to sale by or on the order of a physician. See Instructions for Use for full prescribing information, including indications, contraindications, warnings, precautions, and adverse events.

NK NetworkCare™ is a registered trademark of Nihon Kohden America, LLC.

**Comprehensive maintenance and protection
to ensure uninterrupted patient care.**

NK NetworkCare™ Protection Features

Real-Time Monitoring

- Continuous hardware failure detection
- Error log analysis
- Alarm message monitoring
- Network performance optimization

Security Protection

- Comprehensive malware defense
- Ransomware protection
- Antivirus capabilities
- Breach prevention

System Maintenance

- Network optimization
 - Performance monitoring
 - Hardware maintenance
 - Continuous system updates
-

NK NetworkCare™ Composition

- SolarWinds monitoring platform
- Sumo Logic analytics
- Integrated protection tools
- Real-time monitoring dashboard

**Contact your NK representative to learn how NK Network Care™ can protect and maintain
your facility's critical patient monitoring systems.**

NIHON KOHDEN AMERICA
15353 Barranca Pkwy, Irvine, CA 92618
800.325.0283 us.nihonkohden.com

NK NetworkCare™ is a registered trademark of Nihon Kohden America, LLC.